

References

Honkanen PB, Kellomäki M, Lehtimäki MY et al: Bioreconstructive joint scaffold implant arthroplasty in metacarpophalangeal joints: short-term results of a new treatment concept in rheumatoid arthritis patients. *Tissue Eng.* 2003, 9: 957–65

Honkanen PB, Kellomäki M, Konttinen Y, Mäkelä S, Lehto MA: Midterm Follow-Up Study of Bioreconstructive Poly(lactide Scaffold Implants in Metacarpophalangeal Joint Arthroplasty in Rheumatoid Arthritis Patients. *The Journal of Hand Surgery (European Volume, 2009)* 34E: 2: 179–185

Honkanen, PB, Tiihonen R, Skyttä ET, Ikävalko M, Lehto M, Konttinen Y: Bioreconstructive Poly-L/D-Lactide Implant Compared with Swanson Prosthesis in Metacarpophalangeal Joint Arthroplasty in Rheumatoid Patients: A Randomised Clinical Trial. *J Hand Surg Eur Vol Online First*, published on July 13, 2010 as doi:10.1177/1753193410375777.

Honkanen P: Metacarpophalangeal arthroplasty and partial wrist fusion as a surgical treatment in rheumatoid hand disease. Doctoral thesis. *Acta Universitatis Tamperensis*; 1698; 2012

Ikävalko M, Skyttä E, Belt E: One-Year Results of Use of Poly-L/D-Lactic Acid Joint Scaffolds And Bone Packing in Revision Metacarpophalangeal Arthroplasty. *Journal of Hand Surgery (European Volume, 2007)* 32E: 4: 427–433

Mai S, Mai B: Ein- bis Zwei-Jahreserfahrungen mit einem neuen biodegradierbaren Implantat für kleine Gelenke. *Orthopädische Praxis* 43, 4, 2007: 159-167

Mai B, Mai S: Interposition arthroplasty with a new bioreplaceable implant for small joints RegJoint – five year experience in hand surgery. Abstract FESSH Congress, Antalya, Turkey 2013

Mai S, Hassel K, Mai B: Five year results with a new bioreplaceable implant for MTP joints (RegJoint). Poster EFORT Congress. Istanbul, Turkey 2013.

Mai S, Mai B: Die rheumaorthopädische Behandlung der Hand (Orthopedic Treatment of Rheumatoid Deformed Hand. *OUP* 2013; 9: 400–404. DOI 10.3238/oup.2013.0400–0404

Tiihonen R, Skyttä E, Ikävalko M, Kaarela K, Belt E: Comparison of Bioreplaceable Interposition Arthroplasty with Metatarsal Head Resection in the Rheumatoid Foot. One-Year Results of a Prospective, Randomised Clinical Study. *Foot Ankle Int.* 2010 Jun; 31(6):505-10

Tiihonen R, Skyttä E, Kaarela K, Ikävalko M, Belt E: Reconstruction of the trapeziometacarpal joint in inflammatory joint disease using interposition of autologous tendon or poly-L-D-lactic acid implants: a prospective clinical trial. *J Plast Surg Hand Surg.* 2012 Apr; 46(2):113-9

Tiihonen R: Bioabsorbable Poly-L/D-lactide 96/4 (PLDLA) Implant in Hand and Forefoot Joint in Chronic Inflammatory Arthritis. Doctoral thesis. *Acta Universitatis Tamperensis* 1805. 2013

Complete reference list and copies of references are available from orders@scaffdex.com

B0001RG Rev 7 / 30.4.2016

RegJoint™

- the bioreplaceable small joint spacer

RegJoint™

Mobile life with RegJoint!


RegJoint™ is available in the entire EU area and numerous other countries worldwide. Further information on the product and your local distributor is available on www.scaffdex.com and orders@scaffdex.com.

 scaffdex

Scaffdex Oy
Näsilinnankatu 48 E
FI-33200 Tampere, FINLAND

e-mail: orders@scaffdex.com
www.scaffdex.com

CE 0297





RegJoint™ is a tissue engineered solution for small joint arthrosis in hands and feet. It is a spacer that guides the formation of a functional, flexible fibrous tissue to the joint gap created in joint revision. The new connective tissue will reach an organized structure with a strong collagen framework in three months after the implantation. This neo tissue maintains the function of the mobile joint. RegJoint™ is made of well-known and safe poly-96L/4D-lactide which will totally degrade from the joint cavity in 2-3 years the final degradation products being only water and carbon dioxide.

Advantages

- Significantly decreases pain of operated joints
- Improves functionality and appearance of operated joints
- No weight or function restrictions after rehabilitation period
- Can be used in revision operations of failed non-absorbable implants - also in cases with major diaphyseal osteolysis
- No risk of failure caused by implant breakage

Operational technique

- No specific instruments are needed
- The size of RegJoint™ is chosen to completely cover the resected bone ends (cortical bone)
- Metallic non-disposable RegJoint™ sizers can be used to evaluate the ideal RegJoint™ size
- RegJoint™ can be fixed in its place with resorbable sutures or K-wire
- For further information, please see Instructions for Use and Surgical Technique

Clinical experience

- Thousands of operated implants
- Over 10 years' follow-up times
- Already in clinical trials over 500 joints operated
- Excellent clinical efficacy and high patient satisfaction
- Outstanding safety track record

RegJoint™

product selection



REF number	Height (mm)	Diameter (mm)
RG0001	3.6	8
RG0002	4.0	10
RG0003	4.0	12
RG0004	4.5	14
RG0005	4.5	16
RG0006	4.5	18
RG0007	4.5	20



RegJoint™ sizer

product selection



REF number	Height (mm)	Diameter (mm)
RGS001	3.6	8
RGS002	4.0	10
RGS003	4.0	12
RGS004	4.5	14
RGS005	4.5	16
RGS006	4.5	18
RGS007	4.5	20

