

— THREADED TITANIUM SUTURE ANCHOR

MEMOFIX™



DOCUMENTATION



groupe lépine

Depuis 1714

SUTURE ANCHOR

- titanium aluminium vanadium alloy.
- ready-to-use anchor with holder.
- sterilized by Gamma irradiation.
- implant is double vacuum packed.





SUTURE

- Force Fiber™ braided UHMW-PE.
- for anchor MN2 Ø 2 mm:
 - USP 2-0 (0.35 mm), 1 blue-white strand
- for anchor MN3 Ø 3 mm:
 - USP 2-0 (0.35 mm), 1 blue-white strand and 1 white strand
- for anchors MN4 Ø 4.5 mm and MN5 Ø 5.5 mm:
 - USP 2 (0.6 mm), 1 blue-white strand and 1 white strand
- sterile suture is preassembled on anchor.

NEEDLES

- stainless steel.
- round cross-section.
- for anchor MN2 Ø 2 mm:
 - 2 needles Ø 0.63 mm and length 20 mm.
- for anchor MN3 Ø 3 mm:
 - 4 needles Ø 0.63 mm and length 20 mm.
- for anchors MN4 Ø 4.5 mm and MN5 Ø 5.5 mm:
 - 2 needles Ø 0.98 mm and length 26 mm with 2 needles Ø 0.98 mm and length 36 mm
- sterile needles are preloaded.

MEMOFIX™ threaded titanium suture anchors

<i>Designation</i>	<i>Ref</i>
 SCREWED ANCHOR Ø 2 WIRE Ø 0.35	MN2
 SCREWED ANCHOR Ø 3 WIRE Ø 0.35	MN3
 SCREWED ANCHOR Ø 4.5 WIRE Ø 0.6	MN4
 SCREWED ANCHOR Ø 5.5 WIRE Ø 0.6	MN5

Suture on MEMOFIX™ threaded suture anchors

References	Total length (mm)	Threaded length (mm)	Ø Maximum (mm)
MN2	5,6	3,9	2
MN3	8,8	5,8	3
MN4	11,2	8,2	4,5
MN5	14,8	11,8	5,5

lépine Algérie
9 rue philosophe Tabrizi
Les sources
Bir Mourad Rais - ALGER
lepine-algerie@groupe-lepine.com

lépine Ibérica
J.J. Tadeo Murguía, 3 y 5 Bajos
20304 IRÚN (GUIPÚZCOA)
lepine-iberica@groupe-lepine.com

lépine Colombia
Carrera 49 B #93-49
Barrio la Castellana
110221 BOGOTÁ
lepine-colombia@groupe-lepine.com

lépine Italia
Via Cassanese 100,
20090 SEGRATE (MILANO)
lepine-italia@groupe-lepine.com

lépine Deutschland
Ohle Ring 23-25
STADE-WIEPENKATHEN, 21684
lepine-deutschland@groupe-lepine.com

lépine Maroc
79 avenue IBN SINA
10080 RABAT - AGDAL
lepine-maroc@groupe-lepine.com

175 rue Jacquard - CS 50307
69727 Genay Cedex - FRANCE
Tél. +33 (0)4 72 33 02 95
Fax +33 (0)4 72 35 96 50

www.groupe-lepine.com



groupe lépine

Depuis 1714